

## Digit Widget® Therapy Guidelines

Included with the Digit Widget are forms to graph joint ROM. We recommend that the therapist or surgeon document joint ROM on these graphs at each patient visit, and note the number of hours the device is worn and the strength of the bands used on the forms. This information provides an objective way of determining the effect of the device on the joint.

Instruct the patient to wear the Digit Widget as close to full time as possible, allowing time off for bathing, hand washing, etc. We also recommend that the patient remove the band(s) 4-5 times per day for ~10 minutes to perform un-resisted flexion exercises. Because the bands demonstrate a decrease in tension with use, the bands should be replaced daily. Please review our “Patient Information for the Digit Widget” under the downloads section of our website for proper cuff placement and band changing information.

Start the patient with the lightest weight bands. If the joint demonstrates an increase in extension as a function of time and the patient does not have increasing pain, swelling or stiffness in the finger, the patient can continue using just the light bands. If the finger is not gaining extension, the patient can move up to the medium or heavy bands. The gain in extension from week to week varies with the etiology of the contracture, the length of time it has been present and the number of prior surgeries the finger may have had. On your graphs, look for an extension gain of 10-20° per week without a significant change in flexion.

Some patients do not need to use anything stronger than the light bands, and other patients may require 2-3 heavy bands. It is totally dependent on the joint’s response to torque. Extension torque may require up to 8 weeks to reverse a 90° flexion contracture. Some contractures may reverse more quickly, and some may take longer. Graphing the patient’s ROM during treatment makes it easier to see how joint extension is progressing and to adjust the treatment accordingly.

If the MP joint is hyper-extending, use the MP Flexion Strap that is included inside Pack B. Holding the MP joint in some degree of flexion will help concentrate the extension torque at the PIP joint. Follow the directions for use on the Condensed Assembly Guide or Patient Information Card. The MP Flexion Strap can be trimmed to fit.

If the patient’s skin becomes irritated under the Cuff, a cotton stockinet may be worn under the Cuff.

The device should not be removed entirely until the patient has a trial period without the bands on to determine if the PIP joint has a competent extensor mechanism. This is especially important if the underlying cause of the contracture is still not fully defined. Ideally, the device should be left on for at least 3-4 weeks after maximum extension is achieved. Most patients will require post Digit Widget splinting program for several weeks or even months to maintain gains in active joint extension.

A review of the following resources on our website is also recommended:

FAQ’s <https://handbiolab.com/products/digit-widget/faq/>

Surgical technique video (therapy begins at timestamp 13:06) <https://handbiolab.com/videos/digit-widget/>

If the patient needs more supplies, contact Hand Biomechanics Lab.

Customer Care Team

916-923-5073

website: [www.handbiolab.com](http://www.handbiolab.com)

email: [info@handbiolab.com](mailto:info@handbiolab.com)