

Product Technical Information

Effective: December 5, 2014

Re: Latex Content of Agee Digit Widget[®] System, Item No. DWD-232

The elastic bands provided with the Digit Widget System contain natural rubber latex. These bands are for post surgical use; they are packaged in a plastic vial and are applied by medical staff as a part of the non-surgical assembly. Thereafter the patient changes the bands daily.

For patients with a sensitivity to latex, contact Hand Biomechanics Lab, Inc. at 800-522-5778 to request bands that are 'not made with natural rubber latex'. These alternative bands are weakened with exposure to direct light. It is imperative that your patient understands the importance of changing the bands daily.