

Product Technical Information

Effective: August 14, 2009

Re: Nickel Content of Agee TurnKey FCS® System, Item No FCS-400

Skeletal pins contained in the TurnKey FCS System are made of type 316, implant grade, Stainless steel. This alloy is composed of approximately 14% Nickel, among other metals.

For patients who are allergic to nickel, TurnKey skeletal pins may not be suitable. We do not provide, or know of, an alternative pin for patients who are allergic to nickel.

Contact Hand Biomechanics Lab at 800-522-5778 if you need further technical information on the TurnKey FCS System.