

## **Product Technical Information**

Effective: August 14, 2009

Re: Nickel Content of Agee WristJack® System, Item No.CFD-147

Fixation pins contained in the WristJack System are made of type 316, implant grade stainless steel. This alloy is composed of approximately 14% NICKEL, among other metals.

For patients who are allergic to nickel, WristJack System fixation pins may not be suitable. We do not provide, or know of, an alternative pin for patients who are allergic to nickel.

Contact Hand Biomechanics Lab, Inc. at 800-522-5778 if you need further technical information on the WristJack System.